

# **Specs:**

## **Skill Development:**

- Central Vascular Catheterization, using Subclavian, Supraclavicular, and Internal Jugular approaches
- Using ultrasound guidance during catheter placement
- Palpating anatomic landmarks significant to the procedure

## **Features and Benefits:**

- Excellent for ultrasound guided central venous access hands-on training
- Anatomically correct human torso with landmarks
- Differentiates arterial and venous blood to show positive or negative results
- Self sealing veins and skin for multiple cannulations
- Replaceable tissue sets allow for multiple uses
- Arterial pulse and proper landmarks are present to help avoid and detect errors

## **Product Components:**

- CentraLineMan Torso
- Replaceable Tissue (CLMT-40)
- Vessel Filling Port
- Fill Line
- Pulse - Hand Pump
- Venous Pressure Regulator
- Carrying Case
- Arterial Fluid (red)
- Venous Fluid (blue)
- User's Guide

## **Dimensions:**

- Size: 18" h; 27" w; 8.5" d (carrying case)
- Weight: 27 lbs